



1055

Oxygen therapy mask for adults with bag



FEATURES

High concentration oxygen therapy mask for adults. Produced with non-toxic material, complete with bag, adjustable elastic, clip for nose and connecting tube of cm 210.

MAIN FEATURES

Size

Adults

SEE ALSO



1050

Oxygen therapy mask for adults



1045

Oxygen therapy mask for children
with bag



1040

Oxygen therapy mask for children

ACCESSORY OF



1366

Red knapsack "115" with side
pockets and compartments
COMPLETE
