



## 19130-956

"MTX-956" single space support with lateral movement and ramp of mm. 180, equipped with EN1789 certified fixing system art. 956 - EN1789 certified




### FEATURES

Single space stretcher support with lateral movement and ramp mm 180 in the underneath space can be stored a scoop stretcher or a spinal board. On request supplied with UNI EN 1789 certified blocking devices Art. 955 or 956.



Lateral movement: cm 30

### MAIN FEATURES

Side sliding	
Width (cm)	60,5
Length (cm)	214,5
Weight (kg)	66,5
Loading capacity (kg)	400
Height (cm)	19,5

---

IMAGES

---



---

ACCESSORIES

---



**19114**  
130 mm ramp for stretcher  
supports



**19112**  
100 mm ramp for stretcher  
supports



**19123**  
Ribaltina 130 mm con attacchi per  
supporti e rullo



**19115**  
Ribaltina da 180 mm con attacchi  
per supporti e rullo

---