



19130-956

"MTX-956" single space support with lateral movement and ramp of mm. 180, equipped with EN1789 certified fixing system art. 956 - EN1789 certified



FEATURES

Single space stretcher support with lateral movement and ramp mm 180 in the underneath space can be stored a scoop stretcher or a spinal board. On request supplied with UNI EN 1789 certified blocking devices Art. 955 or 956.



Lateral movement: cm 30

MAIN FEATURES

Side sliding	
Width (cm)	60,5
Length (cm)	214,5
Weight (kg)	66,5
Loading capacity (kg)	400
Height (cm)	19,5

IMAGES



ACCESSORIES



19114
130 mm ramp for stretcher supports



19112
100 mm ramp for stretcher supports



19123
Ribaltina 130 mm con attacchi per supporti e rullo



19115
Ribaltina da 180 mm con attacchi per supporti e rullo
