



19136

"MX" fixed single space stretcher support with ramp of mm 180 - EN1789 certified



FEATURES

Fixed single space stretcher support with ramp mm. 180. On request supplied with EN1789 certified blocking devices art. 955 or 956.



MAIN FEATURES

Width (cm)	60,5
Length (cm)	208,5
Weight (kg)	41,5
Loading capacity (kg)	400
Plate height	12,5
Total height	15,5

ACCESSORIES



19114
130 mm ramp for stretcher supports



19112
100 mm ramp for stretcher supports