



## 19136

"MX" fixed single space stretcher support with ramp of mm 180 - EN1789 certified



### FEATURES

Fixed single space stretcher support with ramp mm. 180. On request supplied with EN1789 certified blocking devices art. 955 or 956.



### MAIN FEATURES

Width (cm)	60,5
Length (cm)	208,5
Weight (kg)	41,5
Loading capacity (kg)	400
Plate height	12,5
Total height	15,5

### ACCESSORIES



**19114**  
130 mm ramp for stretcher supports



**19112**  
100 mm ramp for stretcher supports