



## 19136-955

"MX-955" fixed single space stretcher support with ramp of mm. 180, equipped with certified fixing system art. 955 - EN1789 certified



### FEATURES

Fixed single space stretcher support with ramp mm. 180. On request supplied with UNI EN 1789 certified blocking devices Art. 955.



### MAIN FEATURES

Width (cm)	60,5
Length (cm)	208,5
Weight (kg)	50
Loading capacity (kg)	400
Height (cm)	13,5

### ACCESSORIES



**19114**  
130 mm ramp for stretcher supports



**19112**  
100 mm ramp for stretcher supports