



19140

"MTX-2" double space stretcher support with lateral movement and ramp of mm. 350 - EN1789 certified



FEATURES

Double space stretcher support with lateral movement and ramp mm 350 in the underneath space can be stored both a scoop stretcher and a spinal board. On request supplied with UNI EN 1789 certified blocking devices Art. 955 or 956.
Gas spring included.




Lateral movement: cm 30

MAIN FEATURES

Side sliding	<input checked="" type="checkbox"/>
Width (cm)	63
Length (cm)	229
Weight (kg)	76
Loading capacity (kg)	400
Plate height	25
Total height	28

IMAGES



ACCESSORIES



19110
Stainless steel plate mm 1940 x 572
for stretcher supports



787L
Gas spring for ramp



19123
Ribaltina 130 mm con attacchi per
supporti e rullo



19115
Ribaltina da 180 mm con attacchi
per supporti e rullo

STANDARD ACCESSORIES



821C
Gas spring for ramp
