




## 19140-955

"MTX-2-955" double space stretcher support with lateral movement and ramp of mm. 350, equipped with fixing system art. 955 - EN1789 certified



### FEATURES

Double space stretcher support with lateral movement and ramp mm 350 in the underneath space can be stored both a scoop stretcher and a spinal board. On request supplied with UNI EN 1789 certified blocking devices Art. 955 or 956. Gas spring included.



Lateral movement: cm 30

### MAIN FEATURES

Side sliding	
Width (cm)	63
Length (cm)	229
Weight (kg)	84,5
Loading capacity (kg)	400
Height (cm)	27,5

### IMAGES



---

## ACCESSORIES

---



**19110**  
Stainless steel plate mm 1917 x 572  
for stretcher supports



**787L**  
Gas spring for ramp



**19123**  
Ribaltina 130 mm con attacchi per  
supporti e rullo



**19115**  
Ribaltina da 180 mm con attacchi  
per supporti e rullo

---

## STANDARD ACCESSORIES

---



**821G**  
Gas spring for ramp