




# 19146

"MX-2" double fixed space stretcher support with ramp of mm 350 - EN1789 certified



## FEATURES

Double space fixed stretcher support with ramp mm 350. On request supplied with UNI EN 1789 certified blocking devices Art. 955 or 956.  
Gas spring included.



## MAIN FEATURES

Width (cm)	63
Length (cm)	223
Weight (kg)	59,7
Loading capacity (kg)	400
Plate height	20,5
Total height	23,5

---

## ACCESSORIES

---



**19110**

Stainless steel plate mm 1940 x 572  
for stretcher supports



**787L**

Gas spring for ramp

---

## STANDARD ACCESSORIES

---



**821G**

Gas spring for ramp