




19146-955

"MX-2-955" double space fixed stretcher support with ramp of mm. 350, equipped with certified fixing system art. 955 - EN1789 certified



FEATURES

Double space fixed stretcher support with ramp mm 350. On request supplied with UNI EN 1789 certified blocking devices Art. 955.
Gas spring included.



MAIN FEATURES

| | |
|-----------------------|------|
| Width (cm) | 63 |
| Length (cm) | 223 |
| Weight (kg) | 68,2 |
| Loading capacity (kg) | 400 |
| Height (cm) | 21,5 |

ACCESSORIES



19110
Stainless steel plate mm 1917 x 572
for stretcher supports



787L
Gas spring for ramp

STANDARD ACCESSORIES



821G
Gas spring for ramp
