




## 19146-956

"MX-2-956" double fixed space stretcher support with ramp of mm. 350, equipped with certified fixing system art. 956 - EN1789 certified



### FEATURES

Double space fixed stretcher support with ramp mm. 350. On request supplied with UNI EN 1789 certified blocking devices Art. 956.  
Gas spring included.



### MAIN FEATURES

Width (cm)	63
Length (cm)	223
Weight (kg)	68,2
Loading capacity (kg)	400
Height (cm)	21,5

---

## ACCESSORIES

---



**19110**  
Stainless steel plate mm 1917 x 572  
for stretcher supports



**787L**  
Gas spring for ramp

---

## STANDARD ACCESSORIES

---



**821G**  
Gas spring for ramp

---