



## 1944 Nitrile gloves

---

Throwaway nitrile gloves, size L



---

### FEATURES

Medical throwaway non-sterile nitrile gloves, particularly suitable for people allergic to natural latex rubber and sensitive to type-1 allergic reactions (immediate hypersensitivity). Pale blue, odourless, uniform thickness and free of bubbles, excellent touch and agility sensitivity, they increase comfort even during prolonged operations. Anatomically shaped, ambidextrous with reinforced anti-rollup elastic wristband, no latex protein or residual production chemicals, and without lubricating powder. Available in 4 different sizes, packages containing 100 pcs.

---

### MAIN FEATURES

Material

Nitrile

---

### IMAGES



---

## SEE ALSO

---



**Nitrile gloves**  
Throwaway nitrile gloves, size XL



**Nitrile gloves**  
Throwaway nitrile gloves, size M



**Nitrile gloves**  
Throwaway nitrile gloves, size S

---