



1980

Robin scissors



FEATURES

Rescue scissors perfect for medical emergency operators. Produced with stainless steel and with handle covered by plastic material. It is used to break glass and to cut security belts and clothes in case of emergency. Complete with wall support.

MAIN FEATURES

Length (cm)	23,8
Width (cm)	11,6
Weight (kg)	0,5
Material	Stainless steel

SEE ALSO



1982
"Rescue" emergency scissors

ACCESSORIES



1950/N

Emergency scissors - black

ACCESSORY OF



1366

Red knapsack "115" with side
pockets and compartments
COMPLETE
