



1080

Rescue mask



FEATURES

Indispensable for all rescue operators, this mask avoid any contacts between mouths during emergency reanimations. Made of transparent material, it is composed by a soft cushion and a valve with non-return membrane. It is all contained in a handy, plastic, shockproof case. Color of case may change.

MAIN FEATURES

Size	Adults
Material	PVC

SEE ALSO



1787
Keyring with mouth-to-mouth
mask

ACCESSORY OF



1366
Red knapsack "115" with side
pockets and compartments
COMPLETE
