



410 Neofix

Neoprene splint for forearm



FEATURES

Rigid flexible splints produced with neoprene and complete with an inner aluminium sheet that can be adjusted as needed during the rescue operations. Inner aluminium sheet can be removed for a practical sanitation. The special material used for their production allow to keep constant body temperature, to adapt splints to each limb and to clean and disinfect them easily. 5 available sizes: leg, arm, forearm, wrist and ankle/elbow.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	
Length (cm)	40
Material	Neoprene

SEE ALSO



Halley 886
Vacuum splint set with pump and carrying case



417
"NEOFIX" Neoprene splints set with case