



## 410 Neofix

Neoprene splint for forearm



## **FEATURES**

Rigid flexible splints produced with neoprene and complete with an inner aluminium sheet that can be adjusted as needed during the rescue operations. Inner aluminium sheet can be removed for a practical sanification. The special material used for their production allow to keep constant body temperature, to adapt splints to each limb and to clean and disinfect them easily. 5 available sizes: leg, arm, forearm, wrist and ankle/elbow.

## MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	<b>₩</b>
Length (cm)	40
Material	Neoprene

## **SEE ALSO**



**Halley 886**Vacuum splint set with pump and carrying case



"NEOFIX" Neoprene splints set with case