



## 412 Neofix

Neoprene splint for wrist



### FEATURES

Rigid flexible splints produced with neoprene and complete with an inner aluminium sheet that can be adjusted as needed during the rescue operations. Inner aluminium sheet can be removed for a practical sanification. The special material used for their production allow to keep constant body temperature, to adapt splints to each limb and to clean and disinfect them easily. 5 available sizes: leg, arm, forearm, wrist and ankle/elbow.

### MAIN FEATURES

|             |                |
|-------------|----------------|
| Size        | Adults         |
| Usage       | Immobilization |
| Trauma      |                |
| Length (cm) | 32             |
| Material    | Neoprene       |

### SEE ALSO



**Halley 886**  
Vacuum splint set with pump and carrying case



**417**  
"NEOFIX" Neoprene splints set with case