



## 412 Neofix

---

Neoprene splint for wrist



---

### FEATURES

Rigid flexible splints produced with neoprene and complete with an inner aluminium sheet that can be adjusted as needed during the rescue operations. Inner aluminium sheet can be removed for a practical sanitification. The special material used for their production allow to keep constant body temperature, to adapt splints to each limb and to clean and disinfect them easily. 5 available sizes: leg, arm, forearm, wrist and ankle/elbow.

---

### MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	
Length (cm)	32
Material	Neoprene

---

### SEE ALSO



**Halley 886**  
Vacuum splint set with pump and carrying case



**417**  
"NEOFIX" Neoprene splints set with case