



417

“NEOFIX” Neoprene splints set with case



FEATURES

Rigid flexible splints produced with neoprene and complete with an inner aluminium sheet that can be adjusted as needed during the rescue operations. Inner aluminium sheet can be removed for a practical sanitification. The special material used for their production allow to keep constant body temperature, to adapt splints to each limb and to clean and disinfect them easily. 5 available sizes: leg, arm, forearm, wrist and ankle/elbow.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Trauma	

SEE ALSO



Halley 886

Vacuum splint set with pump and carrying case



Neofix 412

Neoprene splint for wrist



Neofix 410

Neoprene splint for forearm



Neofix 415

Neoprene splint for ankle/elbow



Neofix 407

Neoprene splint for arm



Neofix 405

Neoprene splint for leg
