




## 6100

---

Long folding stretcher




---

### FEATURES

Stretcher that folds in two, equipped with feet and handles for transport. Aluminium frame and PVC impermeable and tear-resistant bed. Ideal for in-surgery transport, easily positioned. Available in two different dimensions.

---

### MAIN FEATURES

Foldable	
Lying surface	Without Trendelenburg
Length (cm)	167
Width (cm)	50
Height (cm)	13
Folded height (cm)	83
Thickness (cm)	13
Weight (kg)	5,5

---

SEE ALSO

---



**6300/B**  
Vacuum mattress for veterinary use



**Maxima 6200**  
Pediatric anti trauma scoop  
stretcher

---