




6150


Short folding stretcher for veterinary use



FEATURES

Stretcher that folds in two, equipped with feet and handles, aluminium frame and PVC impermeable and tear-resistant bed. Ideal for in- surgery transport, easily positioned dimensions.

MAIN FEATURES

Foldable	
Lying surface	Without Trendelenburg
Length (cm)	127
Width (cm)	50
Height (cm)	13
Folded height (cm)	63
Thickness (cm)	13
Weight (kg)	5,5

SEE ALSO



6300/B
Vacuum mattress for veterinary use



Maxima 6200
Pediatric anti trauma scoop
stretcher



6100
Long folding stretcher
