



## 616

---

Certified vertical fastening system for scoop stretchers



---

### FEATURES

Aluminium support with safe quick locking device for wall mount or vertical storage of a scoop stretcher. The stretcher can be folded and then locked.

With our item 616/A it is possible to fix the scoop stretcher folded.



Device certified according to the European harmonized standards UNI EN 1789

---

### MAIN FEATURES

Length (cm)	44
Width (cm)	2,5
Height (cm)	7,5
Weight (kg)	0,5

---

SEE ALSO

---



**Maxima 630**  
Grey anodized anti trauma scoop  
stretcher



**616/A**  
Certified vertical fixing system for  
scoop stretchers



**Maxima 631**  
Grey anodized anti trauma scoop  
stretcher with head/trunk in 1 pc



**Maxima 630/V**  
Green painted anti trauma scoop  
stretcher



**Maxima 630/G**  
Yellow painted anti trauma scoop  
stretcher



**Maxima 630/B**  
Blue painted anti trauma scoop  
stretcher



**Maxima 630/A**  
Orange painted anti trauma scoop  
stretcher

---