



616/A

Certified vertical fixing system for scoop stretchers



FEATURES

Aluminium supports with safe and rapid hook to fix the scoop stretcher in vertical position to the wall or inside the ambulance. With our item 616/A it is possible to fix the scoop stretcher folded.



Device certified according to the European harmonized safety standards UNI EN 1789.

MAIN FEATURES

Length (cm)	44
Width (cm)	2,5
Height (cm)	9,5
Weight (kg)	0,55

SEE ALSO



616

Certified vertical fastening system
for scoop stretchers

ACCESSORIES



Maxima 630/V
Green painted anti trauma scoop
stretcher



Maxima 630/G
Yellow painted anti trauma scoop
stretcher



Maxima 630/B
Blue painted anti trauma scoop
stretcher



Maxima 631
Grey anodized anti trauma scoop
stretcher with head/trunk in 1 pc



Maxima 630
Grey anodized anti trauma scoop
stretcher



Maxima 630/A
Orange painted anti trauma scoop
stretcher
