



6200 Maxima

Pediatric anti trauma scoop stretcher



FEATURES

Completely constructed in aluminium, the structure is a break-apart design, dividing into two parts, allowing the stretcher to be positioned under the patient without having to remove it. It is equipped with two quick-release buckles.

MAIN FEATURES

Size	Children
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	130
Weight (kg)	7,8
Loading capacity (kg)	170
Pediatric	✓
Material	Anodized aluminium

SEE ALSO



Maxima 6200 Pediatric anti trauma scoop stretcher



Maxima 6200/V Scoop stretcher for veterinary use