



## 6200 Maxima

Pediatric anti trauma scoop stretcher



### FEATURES

Completely constructed in aluminium, the structure is a break-apart design, dividing into two parts, allowing the stretcher to be positioned under the patient without having to remove it. It is equipped with two quick-release buckles.

### MAIN FEATURES

Size	Children
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	130
Weight (kg)	7,8
Loading capacity (kg)	170
Pediatric	
Material	Anodized aluminium

---

SEE ALSO

---



**Maxima 6200**  
Pediatric anti trauma scoop  
stretcher



**Maxima 6200/V**  
Scoop stretcher for veterinary use

---