



624 Double


Head immobilizer for scoop stretcher



FEATURES

This head immobilizer can be mounted and separated on the scoop stretcher. It is supplied with the standard side rigid blocks. Supplied with padded belts for the fixing. It is produced with expanded material at closed cells with an inside aluminium sheet and it is covered by a liquid-proof and bacterial-proof material.

MAIN FEATURES

Size	Adults
Usage	Immobilization
Width (cm)	40
Depth (cm.)	25
Height (cm)	15
Weight (kg)	0,7
Trauma	

IMAGES



ACCESSORY OF



Maxima 631

Grey anodized anti trauma scoop
stretcher with head/trunk in 1 pc



Maxima 630

Grey anodized anti trauma scoop
stretcher
