



630 Maxima

Grey anodized anti trauma scoop stretcher



FEATURES

Thanks to its features the scoop stretcher can be placed under the patient with ease, without lifting the body. Built with carefully designed materials, MAXIMA ensures solidity, reliability and durability.

- Efficient splitting system with locking device
- Adjustable length with automatic locking device
- Bag with n. 3 straps (quick release) for the immobilization of the patient (603/MEB)
- Easily arranged inside the ambulance
- Built with light and sturdy aluminium
- Carrying case available on request (Art. 630/CS)



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Folded length (cm)	120
Epaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
10G Certification	<input checked="" type="checkbox"/>
Material	Anodized aluminium
Trauma	<input checked="" type="checkbox"/>
Alpine	<input checked="" type="checkbox"/>

IMAGES



SEE ALSO



Maxima 630/N
Black painted anti trauma scoop
stretcher



Maxima 630/V
Green painted anti trauma scoop
stretcher



Maxima 630/G
Yellow painted anti trauma scoop
stretcher



Maxima 630/B
Blue painted anti trauma scoop
stretcher



Maxima 630/A
Orange painted anti trauma scoop
stretcher



Maxima 631
Grey anodized anti trauma scoop
stretcher with head/trunk in 1 pc

ACCESSORIES



614
Certified horizontal fastening
system for scoop stretchers



616/A
Certified vertical fixing system for
scoop stretchers



616
Certified vertical fastening system
for scoop stretchers



Double 624
Head immobilizer for scoop
stretcher



630/CS
Carrying case for scoop stretchers



603/MEB
Nylon bag containing 3 belts Art.
6077/Ameb

ACCESSORY OF



603/MEB

Nylon bag containing 3 belts Art.
607/Ameb
