



# 630/B Maxima

---

Blue painted anti trauma scoop stretcher



---

## FEATURES

Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.




Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

- Efficient splitting system with locking device
- Adjustable length with automatic locking device
- Bag with n. 3 straps (quick release) for the immobilization of the patient (603/MEB)
- Easily arranged inside the ambulance
- Built with light and sturdy aluminium
- Carrying case available on request (Art. 630/CS)

---

## MAIN FEATURES

---

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Epaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
10G Certification	
Material	Painted aluminium

---

## SEE ALSO

---



### Maxima 631

Grey anodized anti trauma scoop  
stretcher with head/trunk in 1 pc



### Maxima 630/N

Black painted anti trauma scoop  
stretcher



### Maxima 630

Grey anodized anti trauma scoop  
stretcher



### Maxima 630/A

Orange painted anti trauma scoop  
stretcher

---

## ACCESSORIES

---



### 616/A

Certified vertical fixing system for  
scoop stretchers



### 630/CS

Carrying case for scoop stretchers



### 614

Certified horizontal fastening  
system for scoop stretchers



### 616

Certified vertical fastening system  
for scoop stretchers



### Double 624

Head immobilizer for scoop  
stretcher



### 603/MEB

Nylon bag containing 3 belts Art.  
607/Ameb