



# 630/B Maxima

Blue painted anti trauma scoop stretcher



## **FEATURES**

Keeping the very same features of Art.630, thanks to a special paint (designed for optimal durability), this stretcher can be used to differentiate operators or a specific area of operation.



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

- Efficient splitting system with locking device
- Adjustable length with automatic locking device
- Bag with n. 3 straps (quick release) for the immobilization of the patient (603/MEB)
- Easily arranged inside the ambulance
- Built with light and sturdy aluminium
- Carrying case available on request (Art. 630/CS)

## MAIN FEATURES

Size	Adults
Usage	Rescue
Height (cm)	7
Width (cm)	42
Length (cm)	167
Max length (cm)	202
Epaisseur pliée (cm)	9
Weight (kg)	10
Loading capacity (kg)	170
10G Certification	<u> </u>
Material	Painted aluminium

### **SEE ALSO**



**Maxima 631** Grey anodized anti trauma scoop stretcher with head/trunk in 1 pc



**Maxima 630/N** Black painted anti trauma scoop stretcher



**Maxima 630**Grey anodized anti trauma scoop stretcher



**Maxima 630/A**Orange painted anti trauma scoop stretcher

## **ACCESSORIES**





**630/CS**Carrying case for scoop stretchers



Certified horizontal fastening system for scoop stretchers



**616**Certified vertical fastening system for scoop stretchers



**Double 624** Head immobilizer for scoop stretcher



**603/MEB**Nylon bag containing 3 belts Art.
607/Ameb