



## 638 Rescue stretcher

Collapsible rescue stretcher with adjustable backrest



### FEATURES

Thanks to its particular features this stretcher is suitable for auxiliary uses as secondary stretcher inside the ambulance.

- Safe, adjustable backrest operated by a gas spring and lockable in any position
- Frame in light and sturdy aluminium
- Bag with n. 2 quick unhooking belts for the immobilization of the patient (art. 604/MEB)
- It can be folded for an easy arrangement
- Wheels and feet turned inwards in order to reduce the clearance
- Orange nylon sheet fire proof, water proof, anti-tear and easy to clean

### MAIN FEATURES

Lying surface	Without Trendelenburg
Size	Adults
Usage	Transport
Foldable	
Height (cm)	14
Width (cm)	49
Length (cm)	187
Folded length (cm)	98
Weight (kg)	9,2
Loading capacity (kg)	170
Epaisseur pliée (cm)	17
Material	Anodized aluminium Nylon

---

## SEE ALSO

---



**Rescue stretcher 640**  
Collapsible rescue stretcher

---

---

## ACCESSORIES

---



**604/MEB**  
Orange bag with 2 straps with plastic buckle and ME.BER. ribbon



**964**  
Front fastening system for rescue stretcher art. 638 - 639 - 640



**Rear fixing system**  
Rear fastening system for ambulance cots H=129 mm

---

---

## STANDARD ACCESSORIES

---



**604/MEB**  
Orange bag with 2 straps with plastic buckle and ME.BER. ribbon

---