



638 Rescue stretcher

Collapsible rescue stretcher with adjustable backrest



FEATURES

Thanks to its particular features this stretcher is suitable for auxiliary uses as secondary stretcher inside the ambulance.

- Safe, adjustable backrest operated by a gas spring and lockable in any position
- Frame in light and sturdy aluminium
- Bag with n. 2 quick unhooking belts for the immobilization of the patient (art. 604/MEB)
- It can be folded for an easy arrangement
- Wheels and feet turned inwards in order to reduce the clearance
- Orange nylon sheet fire proof, water proof, anti-tear and easy to clean

MAIN FEATURES

Lying surface	Without Trendelenburg
Size	Adults
Usage	Transport
Foldable	
Height (cm)	14
Width (cm)	49
Length (cm)	187
Folded length (cm)	98
Weight (kg)	9,2
Loading capacity (kg)	170
Epaisseur pliée (cm)	17
Material	Anodized aluminium Nylon

SEE ALSO



Rescue stretcher 640
Collapsible rescue stretcher

ACCESSORIES



604/MEB
Orange bag with 2 straps with plastic buckle and ME.BER. ribbon



964
Front fastening system for rescue stretcher art. 638 - 639 - 640



Rear fixing system
Rear fastening system for ambulance cots H=129 mm

STANDARD ACCESSORIES



604/MEB
Orange bag with 2 straps with plastic buckle and ME.BER. ribbon
