




# 640 Rescue stretcher

Collapsible rescue stretcher



### **FEATURES**

Suitable as a secondary stretcher inside of an ambulance. No backrest. Supplied with a bag containing two straps for patient immobilization (Art. 604/MEB).

#### MAIN FEATURES

Lying surface	Without Trendelenburg
Size	Adults
Usage	Transport
Foldable	<b>₩</b>
Height (cm)	14
Width (cm)	48
Length (cm)	186
Epaisseur pliée (cm)	17
Weight (kg)	6,5
Loading capacity (kg)	170
Folded length (cm)	92
Material	Anodized aluminium Nylon

# SEE ALSO



Rescue stretcher 638
Collapsible rescue stretcher with adjustable backrest

## **ACCESSORIES**



964 Front fastening system for rescue stretcher art. 638 - 639 - 640



Rear fixing system
Rear fastening system for ambulance cots H=129 mm

### STANDARD ACCESSORIES



**604/MEB**Orange bag with 2 straps with plastic buckle and ME.BER. ribbon