




640 Rescue stretcher

Collapsible rescue stretcher



FEATURES

Suitable as a secondary stretcher inside of an ambulance. No backrest.
Supplied with a bag containing two straps for patient immobilization (Art. 604/MEB).

MAIN FEATURES

Lying surface	Without Trendelenburg
Size	Adults
Usage	Transport
Foldable	
Height (cm)	14
Width (cm)	48
Length (cm)	186
Epaisseur pliée (cm)	17
Weight (kg)	6,5
Loading capacity (kg)	170
Folded length (cm)	92
Material	Anodized aluminium Nylon

SEE ALSO



Rescue stretcher 638
Collapsible rescue stretcher with
adjustable backrest

ACCESSORIES



964

Front fastening system for rescue
stretcher art. 638 - 639 - 640



Rear fixing system

Rear fastening system for
ambulance cots H=129 mm

STANDARD ACCESSORIES



604/MEB

Orange bag with 2 straps with
plastic buckle and ME.BER. ribbon