



643 Folding stretcher

Collapsible rescue stretcher foldable in 4 pieces



FEATURES

Built in tubular anodized aluminium, canvas in fire-resistant, impermeable, tear-resistant PVC-coated nylon. Fixed handles. Ideal for storage inside of any emergency vehicle. Straps are not included. (see "Restraint systems" Art. 604/MEB).

MAIN FEATURES

Size	Adults
Usage	Transport
Foldable	
Min height (cm)	12
Width (cm)	49
Max length (cm)	195
Folded length (cm)	97
Epaisseur pliée (cm)	17
Weight (kg)	5
Loading capacity (kg)	170
Material	Anodized aluminium Nylon

IMAGES



SEE ALSO



Folding stretcher 646
Collapsible rescue stretcher
foldable in 2 pieces

ACCESSORIES



604/MEB
Orange bag with 2 straps with
plastic buckle and ME.BER. ribbon
