




## 644 Fixed stretcher

---

Stiff fixed stretcher for storage



---

### FEATURES

Fixed rigid stretcher for storage, easily piled on top of same-type stretchers for storage optimization. Supplied with 2 quick-release straps for patient immobilization (art. 604/MEB).

- Light and handy
- Canvas in orange nylon, fire proof, impermeable, tear-resistant and easy to clean
- Sturdy aluminium frame
- Non-foldable

---

### MAIN FEATURES

Lying surface	Without Trendelenburg
Size	Adults
Usage	Transport
Height (cm)	17
Width (cm)	48
Max length (cm)	186
Weight (kg)	4,9
Loading capacity (kg)	170
Material	Anodized aluminium Nylon

---

### IMAGES



---

## SEE ALSO

---



**Rigid stretcher 644/S**  
Stiff fixed stretcher for storage with  
backrest

---

---

## ACCESSORIES

---



**604/MEB**  
Orange bag with 2 straps with  
plastic buckle and ME.BER. ribbon

---

---

## STANDARD ACCESSORIES

---



**604/MEB**  
Orange bag with 2 straps with  
plastic buckle and ME.BER. ribbon

---