



646 Folding stretcher

Collapsible rescue stretcher foldable in 2 pieces



FEATURES

Designed for the transport also of heavy patients, this stretcher can be folded in two pieces. The canvas is built in fire-resistant, impermeable, tear-resistant PVC-coated nylon. Fixed handles. Straps are not included. (see "Restraint systems" Art. 604/MEB).

MAIN FEATURES

Size	Adults
Usage	Transport
Foldable	
Min height (cm)	14
Width (cm)	50
Max length (cm)	201
Folded length (cm)	101
Epaisseur pliée (cm)	17
Weight (kg)	6,2
Loading capacity (kg)	170
Material	Anodized aluminium Nylon

IMAGES



SEE ALSO



Folding stretcher 643
Collapsible rescue stretcher
foldable in 4 pieces

ACCESSORIES



604/MEB
Orange bag with 2 straps with
plastic buckle and ME.BER. ribbon
