



## 650 Maya

Folding sedan chair



### FEATURES

Designed to be an essential device for patient transport this chair is made of light anodized aluminium. Supplied with 4 levers (2 rear and 2 front extractable) and 2 fixed wheels Ø 125 mm, it is practical and easy to use. Open/close mechanism avoiding accidental closures. Supplied complete with black bag with 3 belts (Art. 590/N).



Device certified according to the European harmonized standards UNI EN1865 - UNI EN1789

### MAIN FEATURES

Foldable	
Wheel type	2 Small Fixed Wheels
Type	
Height (cm)	91
Width (cm)	52
Length (cm)	74
Max length (cm)	128
Epaisseur pliée (cm)	24
Weight (kg)	10
Loading capacity (kg)	170
Folded height (cm)	90
10G Certification	
Cardio	
Material	Anodized aluminium

---

## IMAGES

---



---

## SEE ALSO

---



### **Maya 653**

Folding sedan chair with 4 wheels,  
footboards and arms

---

## ACCESSORIES

---



### **618**

Fixing system for sedan chairs

---

---

## STANDARD ACCESSORIES

---



**590/N**

Black bag with 3 belts for sedan  
chairs

---