



650 Maya

Folding sedan chair



FEATURES

Designed to be an essential device for patient transport this chair is made of light anodized aluminium. Supplied with 4 levers (2 rear and 2 front extractable) and 2 fixed wheels Ø 125 mm, it is practical and easy to use. Open/close mechanism avoiding accidental closures. Supplied complete with black bag with 3 belts (Art. 590/N).



Device certified according to the European harmonized standards UNI EN1865 - UNI EN1789

MAIN FEATURES

10G Certification	
Foldable	
Wheel type	2 Small Fixed Wheels
Type	
Height (cm)	91
Width (cm)	52
Length (cm)	74
Max length (cm)	128
Epaisseur pliée (cm)	24
Weight (kg)	10
Loading capacity (kg)	170
Folded height (cm)	90
Cardio	
Seat lenght (cm)	43
Material	Anodized aluminium

IMAGES



SEE ALSO



Maya 653

Folding sedan chair with 4 wheels,
footboards and arms

ACCESSORIES



618

Fixing system for sedan chairs

STANDARD ACCESSORIES



590/N

Black bag with 3 belts for sedan
chairs
