



653 Maya

Folding sedan chair with 4 wheels, footboards and arms



FEATURES

Designed to help rescue personnel during rescue procedures. Made of anodized aluminium. Supplied with 4 wheels (2 rotating with brake \emptyset 100 mm and 2 fixed \emptyset 125 mm), easy and safe open/close system, 4 rear handles and 2 front extractable ones, back rest and seat made of thermoplastic material easy to clean and footrest. Supplied complete with black bag containing 3 belts (Art. 590/N).



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

10G Certification	₩
Foldable	⊌
Footboards and arms	₩
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Type	
Height (cm)	90
Width (cm)	52
Length (cm)	75
Max length (cm)	131
Epaisseur pliée (cm)	24
Weight (kg)	14
Loading capacity (kg)	170
Folded height (cm)	98
Material	Anodized aluminium
Seat lenght (cm)	43
Cardio	₩

IMAGES



SEE ALSO



ACCESSORIES



STANDARD ACCESSORIES



590/N Black bag with 3 belts for sedan chairs