



658/Y Quick

Folding sedan chair with 2 wheels – yellow/black



FEATURES

Designed to safely carry patient, and to be opened/closed in simple and safe way. A special mechanism avoids accidental closure of the chair. Supplied with 2 wheels Ø 125 mm, 2 rear levers and 2 front extractable ones. Easy to use. Back rest and seat made of thermoplastic material easy to clean. Supplied complete with black bag containing 3 belts (Art. 590/N).



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789.

MAIN FEATURES

| | |
|-----------------------|----------------------|
| Foldable | |
| Wheel type | 2 Small Fixed Wheels |
| Type | |
| Height (cm) | 90 |
| Width (cm) | 51 |
| Length (cm) | 71 |
| Max length (cm) | 125 |
| Épaisseur pliée (cm) | 23 |
| Weight (kg) | 10 |
| Loading capacity (kg) | 170 |
| Folded height (cm) | 90 |
| 10G Certification | |
| Cardio | |
| Material | Painted steel |

SEE ALSO



Quick 659/Y

Folding sedan chair with 4 wheels,
footboards and arms - yellow/black

ACCESSORIES



618

Fixing system for sedan chairs



612/AN-2/MEB

Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



612/N-2/MEB

Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs



590/N

Black bag with 3 belts for sedan
chairs

STANDARD ACCESSORIES



590/N

Black bag with 3 belts for sedan
chairs
