



658/Y Quick

Folding sedan chair with 2 wheels – yellow/black



FEATURES

Designed to safely carry patient, and to be opened/closed in simple and safe way. A special mechanism avoids accidental closure of the chair. Supplied with 2 wheels Ø 125 mm, 2 rear levers and 2 front extractable ones. Easy to use. Back rest and seat made of thermoplastic material easy to clean. Supplied complete with black bag containing 3 belts (Art. 590/N).



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789.

MAIN FEATURES

Foldable	
Wheel type	2 Small Fixed Wheels
Type	
Height (cm)	90
Width (cm)	51
Length (cm)	71
Max length (cm)	125
Épaisseur pliée (cm)	23
Weight (kg)	10
Loading capacity (kg)	170
Folded height (cm)	90
10G Certification	
Cardio	
Material	Painted steel

SEE ALSO



Quick 659/Y

Folding sedan chair with 4 wheels,
footboards and arms - yellow/black

ACCESSORIES



618

Fixing system for sedan chairs



612/AN-2/MEB

Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



612/N-2/MEB

Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs



590/N

Black bag with 3 belts for sedan
chairs

STANDARD ACCESSORIES



590/N

Black bag with 3 belts for sedan
chairs
