



659/Y Quick

Folding sedan chair with 4 wheels, footboards and arms - yellow/black



FEATURES

Indispensable device for rescue operator. Supplied with 4 wheels (2 rear fixed Ø 125 mm and 2 front rotating with brake Ø 100 mm), foldable arms and footrest, 4 rear handles and 2 front extractables ones. This chair grants safety and comfort to the patient. Back rest and seat made with thermoplastic material easy to clean. Supplied complete with black bag containing 3 belts (Art. 590/N).



Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789

MAIN FEATURES

Foldable	
Footboards and arms	
Wheel type	2 Rotating Wheels + 2 Small Fixed Wheels
Type	
Height (cm)	90
Width (cm)	51
Length (cm)	71
Max length (cm)	126
Epaisseur pliée (cm)	23
Weight (kg)	14
Loading capacity (kg)	170
Folded height (cm)	99
10G Certification	
Cardio	
Material	Painted steel

SEE ALSO



Quick 658/Y
Folding sedan chair with 2 wheels –
yellow/black

ACCESSORIES



618
Fixing system for sedan chairs



612/AN-2/MEB
Adjustable thoracic abdominal belt
orange/black with MeBer band and
4 points fastening for self loading
stretchers and sedan chairs



612/N-2/MEB
Adjustable thoracic abdominal belt
black with MeBer band with 4
points fastening for self loading
stretchers and sedan chairs



590/N
Black bag with 3 belts for sedan
chairs

STANDARD ACCESSORIES



590/N
Black bag with 3 belts for sedan
chairs
