



660 Folding stretcher

Collapsible stretcher foldable in 4 pieces, telescopic handles



FEATURES

This 4 pieces foldable stretcher is the most complete and versatile in the range. The sturdy tubular aluminium frame makes it light and robust. The stretcher features 4 telescopic handles that allow for a larger support surface. The canvas is in fire-resistant impermeable, tear-resistant, easily cleanable grey nylon. Supplied with 2 quick-release straps for patient immobilization (Art. 605/Nmeb).

MAIN FEATURES

Size	Adults
Usage	Transport
Foldable	
Lying surface	Without Trendelenburg
Min height (cm)	18
Epaisseur pliée (cm)	21
Width (cm)	57
Length (cm)	192
Max length (cm)	216
Folded length (cm)	94
Weight (kg)	9,3
Loading capacity (kg)	170
Material	Anodized aluminium Nylon

IMAGES



SEE ALSO



Folding stretcher 1550
Military folding stretcher with
green canvas

ACCESSORIES



660/CS
Carrying case for collapsible
stretcher Art. 660

STANDARD ACCESSORIES



605/Nmeb
50 mm black belt in one piece with
plastic buckle and ME.BER. ribbon
