



7014/A

4 pcs. Skay fire-proof electrowelded mattress



FEATURES

Mattress for self loading stretcher made with Skay and fireproof padding. Totally electrowelded and water proof, can be used on almost all MeBer self loading stretchers (not suitable for stretcher Evox). Easy to clean and sanitize. 3 belts art. 605/C included. Dimensions: cm $186 \times 47 \times 9$.

*All Pictures shown are for illustration purpose only

MAIN FEATURES

Length (cm)	186
Width (cm)	47
Thickness (cm)	9
Material	Skay

STANDARD ACCESSORIES



605/C Orange belt in one piece with plastic buckle mm 50 and ME.BER. ribbon for proof mattresses

ACCESSORY OF



Fixed stretcher 1515 PROOF Fixed stretcher with rigid bed



Mercury Cinque 7075/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF Certified 5 levels self-loading stretcher with 4 swivel wheels -high version



Mercury Lite Cinque 7095/4RG Proof Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels



Mercury Lite Cinque 7097/4RG Proof Certified aluminium 5 levels self-loading stretcher with 4 swivel wheels - high version



Mercury Lite 7082/4RG Proof Certified aluminium variable height stretcher high version with 4 rotating wheels