



7014/A

4 pcs. Skay fire-proof electrowelded mattress



FEATURES

Mattress for self loading stretcher made with Skay and fireproof padding. Totally electrowelded and water proof, **can be used on almost all MeBer self loading stretchers (not suitable for stretcher Evox).**

Easy to clean and sanitize. 3 belts art. 605/C included.

Dimensions: cm 186 x 47 x 9.

**All Pictures shown are for illustration purpose only*

MAIN FEATURES

Length (cm)	186
Width (cm)	47
Thickness (cm)	9
Material	Skay

STANDARD ACCESSORIES



605/C

Orange belt in one piece with plastic buckle mm 50 and ME.BER. ribbon for proof mattresses

ACCESSORY OF



Fixed stretcher 1515 PROOF
Fixed stretcher with rigid bed



Mercury Cinque 7075/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels



Mercury Cinque 7077/4RG PROOF
Certified 5 levels self-loading
stretcher with 4 swivel wheels -
high version



**Mercury Lite Cinque 7095/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels



**Mercury Lite Cinque 7097/4RG
Proof**
Certified aluminium 5 levels self-
loading stretcher with 4 swivel
wheels - high version



Mercury Lite 7082/4RG Proof
Certified aluminium variable height
stretcher high version with 4
rotating wheels