



# 7070/4RG PROOF Mercury

---

Stainless steel variable height stretcher certified with 4 rotating wheels



---

## FEATURES

Mercury is the ideal for the use of the stretcher in tight and narrow places. Operating and easy command, the two frontal wheels can be unlocked or resetted to thier original configuration. Supplied with belts and mattress.



Device certified according to the European harmonized safety standards UNI EN 1865 - UNI EN 1789.

\*All Pictures shown are for illustration purpose only.

---

## MAIN FEATURES

---

Size	Adults
Usage	Transport
10G Certification	
Multilevel	
Lying surface	With Trendelenburg
Wheel type	4 Big Rotating Wheels
Type	Self loading stretcher
Max height (cm)	95
Intermediate height (cm)	72
Max loading height (cm)	66
Min loading height (cm)	64
Folded height (cm)	31
Width (cm)	57
Length (cm)	197
Length with folded trolley (cm)	165
Weight (kg)	50,6
Loading capacity (kg)	250
Material	Stainless steel ABS HDPE (High-density polyethylene)

---

## IMAGES

---



---

## SEE ALSO

---



**Mercury 7072/4RG PROOF**  
Certified stainless steel variable height stretcher high version with 4 rotating wheels



**Mercury Lite 7080/4RG Proof**  
Ambulance cot in aluminium with variable height and four swivel castors



**Mercury Lite 7082/4RG Proof**  
Certified aluminium variable height stretcher high version with 4 rotating wheels



**Mercury Cinque 7075/4RG PROOF**  
Certified 5 levels self-loading stretcher with 4 swivel wheels



**Mercury Cinque 7077/4RG PROOF**  
Certified 5 levels self-loading stretcher with 4 swivel wheels - high version



**Mercury EL 7092 PROOF**  
Certified self-loading stretcher



**Mercury EL 7094 PROOF**  
Certified self-loading stretcher - high version

---

## ACCESSORIES



**938**  
Oximeb oxygen cylinder holder



**939**  
Instrument support for stretchers



**955 PROOF**  
EN 1789 fixing system for Proof roll-in stretchers



**7014/A**  
4 pcs. Sky fire-proof electrowelded mattress



**608/N-MEB**  
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



**12080**  
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



**12106**  
Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



**12080-001**  
Monting kit for bariatric platform on self loading stretcher



**12018N**  
Telescopic foldable phleboclysis stand for self loading stretchers



**12154**  
"BEAR PODS" Bariatric containment system for self-loading stretchers



**7002**  
Electrowelded 4 pcs black mattress for stretchers



**12070/BR**  
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12070/BG**  
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



**12072/BR**  
"COVER OX 5 BR"- professional cylinder holder for 5 lt oxygen bottles



**12072/BG**  
"COVER OX 5 BG"- professional cylinder holder for 5 lt oxygen bottles



**608/A-MEB**  
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**612/N-2/MEB**  
Adjustable thoracic abdominal belt



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt



**707**  
Red 6-handle fireproof cover with MeBer tape for self-loading



**12170**  
Black mattress for right bariatric wing for self-loading stretchers

black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs

orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs

stretcher




**12172**

Black mattress for left bariatric wing for self-loading stretchers



**12084-001**

Cover for bariatric wings set



**12200**

Bariatric wings set for Mercury self-loading stretchers complete with mattresses, belts with reel and assembly kit

---

## STANDARD ACCESSORIES

---



**612/AN-2/MEB**

Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**12025**

Black multi usage object holder for automatic loading stretcher



**608/A-MEB**

2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)