



7210 Frog

"FROG" yellow self loading stretcher with trendelenburg



FEATURES

Made of welded stainless steel powder coated, this is a really strong and reliable device. The Ø 200 mm wheels (2 fixed with brake and 2 rotating with self-centring system) are made of soft material that allows high comfort level to the patient. This model of stretcher has trendelenburg position to 0° from 15° too.

- Rigid platform made of reinforced thermoplastic material.
- Reclining backrest 0° - 80°, operated by a gas piston, can be locked in any position.
- Folding sides.
- Vertical legs protected by plastic buffers.
- Casters are self-centring so that the stretcher can be positioned in the ambulance without them having to be turned manually.
- Front trolley can be lowered to reduce overall dimensions of the stretcher.

MAIN FEATURES

Usage	Transport
Type	Self loading stretcher
Lying surface	With Trendelenburg
Wheel type	2 Rotating Wheels + 2 Big Fixed Wheels
Height (cm)	92
Max loading height (cm)	66
Min loading height (cm)	60
Folded height (cm)	32,5
Width (cm)	56
Length (cm)	195
Length with folded trolley (cm)	169
Loading capacity (kg)	250
Weight (kg)	43
Material	Painted steel ABS HDPE (High-density polyethylene)
Filter feature	

SEE ALSO



Frog 7210 PROOF
"FROG" yellow self loading stretcher
with trendelenburg certified

ACCESSORIES



Front blocking device
Front blocking device for automatic loading stretcher



Rear fixing system
Rear fixing system for roll-in H=148 mm



936
Phleboclysis stand for automatic loading stretcher Winner and Frog



939
Instrument support for stretchers



938
Oximeb oxygen cylinder holder



608/A-MEB
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



608/N-MEB
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



12080
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



12080-001
Monting kit for bariatric platform on self loading stretcher



12025
Black multi usage object holder for automatic loading stretcher



12084
Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



7013/A
3 pcs. Skay fire-proof electrowelded mattress



7013/N
3 pcs. Skay fire-proof electrowelded mattress



12154
"BEAR PODS" Bariatric containment system for self-loading stretchers



12106
Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



12018N
Telescopic foldable phleboclysis stand for self loading stretchers



12070/BG
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



12070/BR
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



12072/BG
"COVER OX 5 BG"- professional
cylinder holder for 5 lt oxygen
bottles



12072/BR
"COVER OX 5 BR"- professional
cylinder holder for 5 lt oxygen
bottles
