



## 7210 PROOF Frog

"FROG" yellow self loading stretcher with trendelenburg certified



### FEATURES

Made of welded stainless steel powder coated, this is a really strong and reliable device. The Ø 200 mm wheels (2 fixed with brake and 2 rotating with self-centring system) are made of soft material that allows high comfort level to the patient. This model of stretcher has trendelenburg position to 0° from 20° too.

- Rigid platform made of reinforced thermoplastic material.
- Reclining backrest 0°- 80°, operated by a gas piston, can be locked in any position.
- Folding sides.
- Vertical legs protected by plastic buffers.
- Casters are self-centring so that the stretcher can be positioned in the ambulance without them having to be turned manually.
- Front trolley can be lowered to reduce overall dimensions of the stretcher.

In coupling with the proper fixation Art. 955 Proof, it is 10g resistant.

This stretcher is supplied complete with mattress (Art. 7013) and belts (nr. 2 Art. 608-A/MEB o 608-N/MEB and nr. 1 Art. 612/AN-2 o Art. 612/N-2).




Device certified according to the European harmonized standards UNI EN 1865 - UNI EN 1789.

---

## MAIN FEATURES

---

Size	Adults
Usage	Transport
10G Certification	
Lying surface	With Trendelenburg
Wheel type	2 Rotating Wheels + 2 Big Fixed Wheels
Type	Self loading stretcher
Height (cm)	92
Max loading height (cm)	66
Min loading height (cm)	60
Folded height (cm)	32,5
Width (cm)	56
Length (cm)	195
Length with folded trolley (cm)	166
Weight (kg)	43
Loading capacity (kg)	250
Material	Painted steel
Filter feature	

---

## SEE ALSO

---



**7265 Proof**  
Aluminium self loading stretcher  
"Frog lite Advance" with  
extractable handles and  
Trendelenburg



**Frog Lite 7260 PROOF**  
"FROG LITE" self loading aluminium  
stretcher with trendelenburg  
certified



**Frog 7210**  
"FROG" yellow self loading stretcher  
with trendelenburg



**Frog 7240 PROOF**  
"FROG PLUS" self loading stretcher  
with trendelenburg certified

---

## ACCESSORIES



**955 PROOF**  
EN 1789 fixing system for Proof roll-in stretchers



**936**  
Phleboclysis stand for automatic loading stretcher Winner and Frog



**939**  
Instrument support for stretchers



**938**  
Oximeb oxygen cylinder holder



**12025**  
Black multi usage object holder for automatic loading stretcher



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt orange/black with MeBer band and 4 points fastening for self loading stretchers and sedan chairs



**612/N-2/MEB**  
Adjustable thoracic abdominal belt black with MeBer band with 4 points fastening for self loading stretchers and sedan chairs



**608/A-MEB**  
2 pieces orange belt with metallic hook and MeBer band (L. 155 cm)



**608/N-MEB**  
2 pieces black belt with metallic hook and MeBer band (L. 155 cm)



**12084**  
Bariatric wings set for self-loading stretchers complete with mattresses and belts with reel



**7013/A**  
3 pcs. Skay fire-proof electrowelded mattress



**7013/N**  
3 pcs. Skay fire-proof electrowelded mattress



**12154**  
"BEAR PODS" Bariatric containment system for self-loading stretchers



**12106**  
Mattress for AIR-COM stretchers - Anti-bedsore mattress with alternate segments



**12018N**  
Telescopic foldable phleboclysis stand for self loading stretchers



**12080**  
"MPB" wide bariatric platform for self-loading stretchers Frog Mercury and Gecko complete of belts with reel and fixing for stretcher



**12080-001**  
Monting kit for bariatric platform on self loading stretcher



**12070/BG**  
"COVER OX 2-3 BG"- professional cylinder holder for 2 and 3 lt oxygen bottles



**12070/BR**  
COVER OX2-3 BR professional cylinder holder with belt and orange MeBer logo for 2 and 3 lt oxygen bottles



**12072/BG**  
"COVER OX 5 BG"- professional  
cylinder holder for 5 lt oxygen  
bottles



**12072/BR**  
"COVER OX 5 BR"- professional  
cylinder holder for 5 lt oxygen  
bottles

## STANDARD ACCESSORIES



**7013/A**  
3 pcs. Skay fire-proof electrowelded  
mattress



**612/AN-2/MEB**  
Adjustable thoracic abdominal belt  
orange/black with MeBer band and  
4 points fastening for self loading  
stretchers and sedan chairs



**608/A-MEB**  
2 pieces orange belt with metallic  
hook and MeBer band (L. 155 cm)