



8004

Colored aneroid sphygmomanometer with carrying case



FEATURES

Range 0-300 mm Hg. Accuracy +/-3 mm Hg. Manometer NON STOP PIN. Adult Velcro® cuff and metal ring. Standard model in vinyl zipper case.

MAIN FEATURES

Size

Adults

SEE ALSO



1740
Rotating wall aneroid
sphygmomanometer



8012
Red palm-style aneroid
sphygmomanometer plus

ACCESSORIES



8143
Double tube pediatric cuff



8142
Single tube pediatric cuff



8140
Adults cuff with single tube



8141
Adults cuff with double tube
