



8004

Colored aneroid sphygmomanometer with carrying case



FEATURES

Sphygmomanometers with nylon cuff with velcro closure, reliable manometer with clip to be fixed on the cuff. Scale 0-300 mmHg. From 260 mmHg to 15 mmHg max deflation time 10 sec. Chromium-plated air tap with air release regulation by screw valve. Supplied with a simulated-leather soft pouch with zip. Latex free.

Range 0-300 mm Hg. Accuracy +/-3 mm Hg. Manometer NON STOP PIN. Adult Velcro® cuff and metal ring. Standard model in vinyl zipper case.

MAIN FEATURES

Siz

Adults

SEE ALSO



1740 Rotating wall aneroid sphygmomanometer



8012Red palm-style aneroid sphygmomanometer plus

ACCESSORIES







8142 Single tube pediatric cuff



